



**ReSkin® Bi-elastic plaster:** Prevents friction  
 | Breathable and kind to your skin | Pain-free removal |  
 Washable and re-useable | Can be trimmed to fit

## ReSkin® Bike Patch

### Protects you while cycling.

Pain due to friction from your cycling shorts while cycling is a phenomena that is very common among cyclists. This can seriously affect the performance of professional cyclists. If this pain can be prevented, performance is bound to increase. Not only that, the pleasure of the amateur biker will increase too.


ReSkin® is the latest innovation from Bioracer that solves this problem. Combining a specially perforated type of Lycra® with a specially created coating has resulted in a skin-friendly product with outstanding qualities. Our product ReSkin® has been put on the market only after extensive testing in both the laboratory and the field. World class cyclists are already using our product with great success. Also very useful for amateurs!

- ReSkin® is applied to the perineum. The perineum is the piece of skin between the anus and the genitals. The perineal area contains sweat glands, hair follicles and sebaceous glands that can often become irritated by rubbing, leading to an infection.
- We advise ladies to cut the strip off along the dotted line marked on the product for enhanced comfort while using ReSkin®.



### Instructions:

1. Remove the backing sheet from the adhesive side of the ReSkin®. Ladies: cut the ReSkin® along the dotted line using scissors. The dotted line is marked on the backing sheet of ReSkin®.
2. Apply the ReSkin® with the sticky side up on top of your finger tips.
3. Bring your hand with the ReSkin® to just underneath the perineum (the piece of skin between the anus and the genitals). Ladies: the ReSkin® can be used both lengthwise and widthwise - whichever you prefer.
4. Place the ReSkin® on the perineum.
5. Rub the remainder of ReSkin® quite firmly onto the inside of your legs.
6. ReSkin® will now protect you from any friction.



### The following applies to all types of ReSkin® products:

- Apply ReSkin® directly to the skin wherever rubbing causes problems. Apply ReSkin® to clean dry skin. **Cut the ReSkin® to the correct size to get a good adhesive bond between the skin and the product.** ReSkin® will immediately relieve pain in areas where rubbed skin is a problem.
- The bi-elasticity of the product (two-directional stretch) creates relief from rubbing by transferring it. Friction now takes place between the ReSkin® and the other materials rather than on the skin. It therefore acts like a second artificial layer of skin.
- ReSkin® is breathable, skin-friendly and painless to remove. The specially created adhesive barely removes a single hair from your body. **Shaving in advance is recommended.**
- The special adhesive used on the surface of the ReSkin® products conforms to ISO 10993-1 and has been tested for cytotoxicity, skin sensitization, and intracutaneous irritation.

**Caution:** ReSkin® is for external use only. Keep out of reach of children. Remove ReSkin® if you experience irritation or discomfort. Please consult your doctor if the irritation or discomfort continues after the product is removed. If you are diabetic or have blood circulation problems please consult your doctor before use.

**Cleaning and instructions for re-use:** Depending on how intensively the ReSkin® is used, it is possible to re-use it several times. When showering after sports, simply stick ReSkin® on the tiles of the bathroom or another suitable smooth area. Rinse the ReSkin® with lukewarm water using the shower head. Use only a mild soap to wash it. Allow the ReSkin® to dry on the surface where it was cleaned. Once it has dried, you can use it again. Tip: Use cling film as a backing sheet when storing it.

International patent pending.  
 ReSkin® is a registered trademark of Bioracer.

INGREDIENTS: Lycra: PES(35%), PUE (10%), Si (55%)  
 Coating: ISO 10993-1, sensitisation intracutaneous irritation tested.

Distributed by/Distributé par: Bescot Healthcare Canada Inc  
 115-639 Dupont St, Toronto, ON, M6G 1Z4

www.ReSkin.ca